



FLYING WITH YOUR

# Tandem Insulin Pump



If you plan to travel by air with your insulin pump, some advance planning may make your security screening go smoother.

Your device is safe for use on commercial airlines and complies with FAA wireless transmission standards.\* It is also designed to withstand common electromagnetic interference and can be safely carried through metal detectors.

However, your Tandem Diabetes Care insulin pump should NOT be put through machines that use X-rays, including airline luggage X-ray machines and full-body scanners.

We recommend disconnecting at the infusion site and asking the security agent for an alternative screening method. If you prefer to stay connected, you can notify the agent about your pump and request to go through a standard metal detector wearing your pump.

It may be helpful to provide a printout of the information below to the security agent during the screening.

If you have any questions about traveling with your insulin pump, please visit [tandemdiabetes.com](https://tandemdiabetes.com) or call Tandem Diabetes Care Customer Support at [\(877\) 801-6901](tel:877-801-6901). We are here for you 24 hours a day, 7 days a week.

\* The pump complies with radiated emissions requirements defined in RTCA/DO-160G, Section 21, Category M.

Fold Here

Show this document to the TSA or airport representative

## Dear TSA or airport representative,

This is an insulin pump from Tandem Diabetes Care, a medically necessary medical device prescribed by a physician. It is designed to withstand common electromagnetic interference and can be safely carried through metal detectors, but it should not be exposed to any form of X-rays. This includes airline luggage X-ray machines and full-body scanners.

Please provide a screening method other than X-ray for this insulin pump. If you have any questions about this device, please call [Tandem Customer Support at \(877\) 801-6901](tel:877-801-6901).


Sincerely,  
Tandem Diabetes Care



**Metal detectors**



**X-rays, including full body scanners and luggage screening machines**



© 2024 Tandem Diabetes Care, Inc. All rights reserved. Tandem Diabetes Care, the Tandem logo, t:slim X2, and Tandem Mobi are either registered trademarks or trademarks of Tandem Diabetes Care, Inc. in the United States and/or other countries. ML-1000524\_E



(877) 801-6901  
[tandemdiabetes.com](https://tandemdiabetes.com)



Need additional help?  
Visit our support center.