



TANDEM MOBI SYSTEM TRAINING Supply Checklist






INSTRUCTIONS

Check off all supplies on this list before your scheduled training.

Tandem Mobi System Supplies



Compatible smartphone with the Tandem Mobi mobile app* installed

Note: Bring phone charging cable.



Tandem Mobi system box



Tandem Mobi system
Accessories & Instructions box




2-3 cartridges



Infusion sets

Note: Bring 2-3 of each infusion set type that you received.



Adhesive sleeves

Continue on next panel →

Additional Supplies



Full vial of Humalog or NovoLog U-100 insulin



Alcohol swabs or skin-prep wipes

Glucose Monitoring Supplies† (all may not apply)



Dexcom G6 CGM

Note: An active sensor can pair to your insulin pump during training.



Blood glucose meter with test strips

* The mobile app requires a compatible iPhone model and operating system (sold separately). Only available to pump users who reside in the United States. A complete list of compatible iPhones will be released closer to commercial availability. † Glucose monitoring supplies sold separately.

Important Safety Information: RX ONLY. Indications for Use: Tandem Mobi system: The Tandem Mobi insulin pump with interoperable technology (the pump) is intended for the subcutaneous delivery of insulin, at set and variable rates, for the management of diabetes mellitus in persons requiring insulin. The pump is able to reliably and securely communicate with compatible, digitally connected devices, including automated insulin dosing software, to receive, execute, and confirm commands from these devices. The pump is intended for single patient, home use and requires a prescription. The pump is indicated for use in individuals 6 years of age and greater.

Warning: Do not use the Tandem Mobi insulin pump before consulting with your healthcare provider. Incorrect settings can result in over-delivery of insulin (which could lead to hypoglycemia) or under-delivery of insulin (which could lead to hyperglycemia). Visit tandemdiabetes.com/safetyinfo for additional important safety information.

© 2024 Tandem Diabetes Care, Inc. All rights reserved. Tandem Diabetes Care, the Tandem logo, Control-IQ, and Tandem Mobi are either registered trademarks or trademarks of Tandem Diabetes Care, Inc. in the United States and/or other countries. Dexcom, Dexcom G6, and any related logos and design marks are either registered trademarks or trademarks of Dexcom, Inc. in the United States and/or other countries. iPhone is a trademark of Apple Inc., registered in the U.S. and other countries. All third-party marks are the property of their respective owners. ML-1013280_A



TANDEM
Diabetes Care

(877) 801-6901
tandemdiabetes.com

FOR USE IN THE
UNITED STATES ONLY

