

Infusion Sets

Finding the right infusion set is an important part of successful insulin pump therapy. That's why we offer so many choices.

AutoSoft 90 Infusion Set





- + 90-degree soft cannula infusion set
- + All-in-one infusion set with preloaded automatic spring inserter makes one-hand insertions quick and easy
- + Large grip area on connector for easier disconnection
- + Reversible connector makes reconnection simple



Preloaded insertion device

Cannula Length	Tubing Length	Cannula Fill	Available Colors
6 mm	23" or 43"	0.3u	Gray, Pink, or Blue
9 mm	23" or 43"	0.5u	Gray, Pink, or Blue

AutoSoft XC Infusion Set




- + 90-degree soft cannula infusion set
- + All-in-one infusion set with preloaded automatic spring inserter makes one-hand insertions quick and easy
- + Reversible connector makes reconnection simple



Preloaded insertion device

Cannula Length	Tubing Length	Cannula Fill	Available Colors
6 mm	23", 32", or 43"	0.3u	Gray, Pink, or Blue
9 mm	23" or 43"	0.5u	Gray

AutoSoft 30 Infusion Set




- + 30-degree soft cannula infusion set
- + All-in-one infusion set with disposable automatic insertion device
- + Reversible connector makes reconnection simple
- + Transparent window at insertion point
- + Often used by active individuals, lean adults, and children



Preloaded insertion device

Cannula Length	Tubing Length	Cannula Fill	Available Colors
13 mm	23" or 43"	0.7u	Gray

VariSoft Infusion Set



- + Allows users to choose their most comfortable insertion angle (20-45 degrees)
- + Easy to insert manually and convenient to detach
- + Transparent window at insertion point
- + Large adhesive pad for extra security against cannula dislodging
- + Often used by active individuals, lean adults, and children

Cannula Length	Tubing Length	Cannula Fill	Available Colors
13 mm	23" or 32"	0.7u	Clear
17 mm	23" or 43"	0.7u	Clear

TruSteel Infusion Set



- + 90-degree stainless steel infusion set
- + Very fine, 29-gauge needle
- + Easy to insert manually and convenient to detach
- + Second adhesive pad for extra security against needle dislodging
- + Simplicity and security for those who have reactions to plastic cannulas and those with a history of bent cannulas
- + Reliability without compromising comfort

Cannula Length	Tubing Length	Cannula Fill	Available Colors
6 mm	23" or 32"	Not Required	Clear
8 mm	23" or 32"	Not Required	Clear



For more information about Tandem's infusion sets and a comparison table of commonly used infusion sets, please visit tandemdiabetes.com/products/infusion-sets

Important Safety Information: RX ONLY. These infusion sets are indicated for the subcutaneous infusion of insulin administered by the t:slim X2 insulin pump for the treatment of diabetes. These infusion sets are indicated for single use. For contraindications, warnings, precautions, and other important information, please refer to the instructions for use accompanying the infusion sets.

© 2022 Tandem Diabetes Care, Inc. All rights reserved. Tandem Diabetes Care, AutoSoft, t:slim X2, TruSteel, and VariSoft are either registered trademarks or trademarks of Tandem Diabetes Care, Inc. in the United States and/or other countries. ML-1001846_E